










# Sterimatic - the safer, cleaner way to inject all Livestock



## Multi-dose system



## Automatic Needle Protector and Cleaning Systems for livestock injections & vaccinations.

It is a widely accepted fact that **if injections are given without changing or cleaning the needle**, there is the potential risk of:

- Bacterial contamination resulting in infection or abscessing and therefore discarded meat.
- Disease transfer between livestock
- Needle-stick or scratch injuries to the operator.

Also, if a dirty needle is used to withdraw a dose from a bottle of injectable, it may contaminate the bottle contents.

### Recommendations:


The Stericap and needle should be replaced after 100 injections.

The Stericap should be used within 3 days from opening.

For more information and data, please visit the Sterimatic website at:

<http://www.sterimatic.com>

## Single-dose system



## Ordering Codes:

Code	Description
STRPA035	Sterimatic 500 dose starter kit, <i>clear</i> , 12mm
STRPA036	Sterimatic 500 dose starter kit, <i>green</i> , 22mm
STRPA037	Sterimatic 500 dose starter kit, <i>blue</i> , 37mm
STRPA001	Refill 500 dose pack for <i>clear</i> (caps and needles)
STRPA002	Refill 500 dose pack for <i>green</i> (caps and needles)
STRPA004	Refill 500 dose pack for <i>blue</i> (caps and needles)
STBAA101	Bottle protector pack for 100ml bottle
STBAA201	Bottle protector pack for 250/500ml bottle

Notes:

Bottle protector packs are 100 dose units